

Fluoroscopy Procedure Table

SONESTA 6210



Features & benefits of the 6210:

Positioning

- Three motors offer unparalleled positioning and access options via a hand held control.
- Four double wheel rolling casters offer unencumbered maneuverability.
- Versatile - eight positioning functions offer the physician standing, seated and supine positioning capabilities.

Imaging

- The side mounted cantilevered design offers the physician anterior / posterior, lateral and oblique views without obstruction.
- Radiolucent carbon fiber seat and back for head-to-toe imaging.
- Provide easier access to, and visualization of, the perineal area.

Functionality

- Allow the patient to remain in one place for preparation, filling, and voiding phases of the urodynamic test.
- Eliminate the risk of catheters falling out or getting misplaced during transfer or ambulation.
- Is rated IPX6 "protection against water incursion".

Procedures commonly performed in the 6210:

- Video Urodynamics
- Imaging Examinations
- Fluoroscopy Procedures

as well as:

- Pelvic Muscle Rehabilitation
- Anorectal Manometry
- Pelvic Exam
- Cystoscopy
- Ultrasound
- Vasectomy
- Colposcopy
- Biopsy
- Colonoscopy
- Gynecologic Procedure



Specifications

Length without foot rest/leg extension	156 cm	63"
Length with foot rest/leg extension	196 cm	78"
Width (total)	99 cm	39"
Height with seat cushion (min.)	50 cm	20"
Height with seat cushion (max.)	130 cm	51"
Imaging area without foot rest/leg extension (length)	156 cm	63"
Imaging area with foot rest/leg extension (length)	196 cm	78"
Weight capacity	250 kg	551 lbs.
Weight capacity, standing position	250 kg	551 lbs.

Electrical

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601
 All products carry the CE mark

Standard Accessories:

1 Hand control	524-6356
1 Plastic funnel with drainage	524-6325-9
1 Arm rest holder	524-6338-4
1 Arm rest	524-6021-6
2 Leg rest holders	525-6338
2 Leg rests	525-6336
2 Fold up foot rests	525-6331-10
1 Side rail	524-6329

