



# Procedure and Examination Chairs

---

## ACCESSORY CATALOG



Model

| 6210 |

Sonesta manufactures and sells tables and chairs for urology, gynecology, urodynamic diagnostics and upper GI procedures.

# Accessories: Sonesta 6210



## Leg & foot rests

6210 Art. number

### Leg rest holder (single)


This leg rest holder is standard on all models. By loosening one handle, you can adjust the leg rest to any height, width or angle. Can be rotated 360 degrees.



525-6338

### Heavy duty rissler leg rests (pair)


This leg rest has a ball and socket construction allowing tilt of the leg cushion in any direction. In addition, it has infinite angling possibilities when mounted in our leg rest holder model 6338. Made of stainless steel.



525-6336

### Fold up footrests (pair)


These footrests are adjustable in length from 30—42 cm/11,8"—16,5". They are also adjustable in width and slide outward 19 cm/7,5" each. The bottom is in the same upholstery as your chair. When folded they function as an extension cushion that transforms your chair into a table where the patient can lie down. The gas piston easily adjusts the angle. They are easily removed with a click-off system.



525-6331

### Footplates (single)

Made of stainless steel with a plastic foot section, fits into our standard leg rest holder.



524-P6413

### Stirrups and clamps (pair)

Stirrups that provide safe and quick positioning which improves surgical site access. The stirrups comes with silicone straps for easier cleaning and attachment. Lithotomy range: +85° to +55°, abduction range: +25° to -9°, boot rotational angle: +/- 25°.



Stirrups & clamps 525-3100

### Benefits & use:

- Multiple positioning options
- Functional design
- Enables optimal examination access
- Engineered ergonomics



Standard configuraton



Optional accessory

# Accessories: Sonesta 6210




## Arm rests

6210

Art. number

### Detachable side rail

Fits on the table top's upper rail free part.  
Used to attach accessories.



S

524-6329

### Arm rest holder

This arm rest holder is used together with 524-6021-6 Arm Rest and 524-6329 detachable side rail. By loosening one handle, you can slide the holder up and down the length of the side rail..




S

524-6338-4

### Arm rest

This arm rest is used together with the 524-6329 detachable side rail and 6338-4 arm rest holder. The arm rest is positioned on the left side of the table and the arm rest on the right side is detachable. The arm rest is covered with a soft foam material.



S

524-6021-6

#### Benefits & use:

- Easy access
- Practical use
- Patient comfort

S

Standard configuration

A

Optional accessory

# Accessories: Sonesta 6210




## Other

6210

Art. number

### Hand control

Waterproof hand control for 6210 table, with cable to transformer.




S

524-6356

### Funnel kit

A removable plastic funnel (top diameter 236 mm/9") with drainage. It is mounted on a fixed height 'quick change' holder.




S

524-6325-9

### Restraint strap

Velcro straps for patient fixation. Pair of 100 cm/39" long and 10 cm/3,9" wide straps.



A

524-P6415



An ADDvise Group Company

HA-2022-05-10  
Version 2

SONESTA MEDICAL AB  
INDUSTRIVÄGEN 7  
171 48 SOLNA, SWEDEN  
TEL: +46 8 502 572 80  
INFO@SONESTAMEDICAL.SE  
ORDER@SONESTAMEDICAL.SE  
WWW.SONESTAMEDICAL.SE

SONESTA MEDICAL INC.  
2 DEBUSH AVENUE UNIT C3  
MIDDLETON, MA 01949, USA  
TEL: +1 630-519-3450  
INFOUS@SONESTAMEDICAL.SE  
ORDERUS@SONESTAMEDICAL.SE  
WWW.SONESTAMEDICAL.SE

S

Standard configuration

A

Optional accessory